

PCR and T-RFLP protocol for BrdU-DNA

*Dr. Dad Roux-Michollet
3516 Bren Hall
University of California Santa Barbara
Santa Barbara, CA 93106-5131*

PCR protocol with f-Taq from Feldan Bio

Master Mix:

Water (qf 50 µl):	23.3 µl	
Buffer (10X):	5 µl	1X (1.5 mM MgCl ₂ and 0.1 mg/ml BSA)
Band Sharpener (5X):	10 µl	1X
MgCl ₂ (25 mM):	1.4 µl	0.7 mM
BSA (10 mg/ml):	1 µl	0.2 mg/ml
dNTP (10 mM each):	1 µl	0.2 mM each
8f-Hex:	2.5 µl	0.5 µM
1389r:	2.5 µl	0.5 µM
Taq (5U/ µl):	0.3 µl	1.5 units

47 µl Master Mix + 3 µl template DNA

Cycling:

- 95°C for 3min
 - 94°C for 30s
 - 53°C for 1min
 - 72°C for 2min
 - 72°C for 10min
- 30 cycles

Supplier:

dNTP (Bioline, Tauton, MA) or Feldan Bio
Primers (Operon Technologies, Alameda, CA),
BSA (New England BioLabs)
MgCl₂ (Bioline)
Taq polymerase (Feldan Bio, Quebec, Canada)

8f: 5'-AGA GTT TGA TCC TGG CTC AG, (Liu et al., 1997)
1389r: 5'-ACG GGC GGT GTG TAC AAG (Osborn et al., 2000)

PCR products purification and quantification

QIAquick PCR purification kit (Qiagen), elution in 50 μ l
Quant-iT dsDNA Assay Kit (Molecular Probes, Invitrogen)

T-RFLP digestion with HhaI (New England BioLabs, Ipswich, MA)

Master Mix:

Water (qf 20 μ l):	X μ l	
NE Buffer 4 (10X):	2 μ l	1X
BSA (100X):	0.2 μ l	1X
HhaI (20 U/ μ l):	0.2 μ l	2 units

Set up reaction in a 96-wells PCR plate

Master Mix: 2.4 μ l to 11 μ l

Purified PCR products: 200 ng (about 17.6 to 9 μ l, depending on the DNA concentration)

Digestion:

37°C for 3h

80°C for 20 min

Cool down at 4°C and quick spin

Freeze at -80°C and send on dry ice to Michigan Sate University

Alternative: Cool down at 4°C and quick spin. Mix 2 μ l of sample digest with 12 μ l of formamide + 0.5 μ l of Map Marker 1000 (standard), then incubate at 94°C for 5 min. It can be stored few days at 4°C before sending to UC Berkeley.

June 17, 2009